

CRYON-X PRO

- Hot/Cold Contrast Therapy
 - Enhanced Compression
 - Rechargeable
 - Portable

The Mobile Solution for Advanced Recovery

→) (←

40

The Cryon-X Pro is a portable therapy system combining cold, heat, and compression in one device, perfect for clinics, athletic facilities, or home use. With a rechargeable battery and software-controlled functionality, it delivers iceless cold therapy and real-time temperature control for targeted pain relief and swelling reduction. Designed for mobility and ease, it's an ideal solution for treating post-surgical and post-traumatic conditions wherever care is needed.



QUICK START GUIDE



Fill device with distilled water as indicated



Connect device to power outlet



Connect one end of hose to device & other to the wrap



Press the power button.



Ensure highlighted steps are completed.



Adjust compression/temperature then click "Start Therapy".



Adjust temperature and save. Then click "START THERAPY."





Adjust compression and save. Then click "START THERAPY."



View website

Cryon-X

25391 Commercentre Dr, Suite 250 Lake Forest, CA 92630

916.587.0295 | info@cryon-x.com | www.cryon-x.com



