



CRYOTHERAPY

Cold Therapy and Pneumatic
Compression for Faster Recovery

POWERFUL RELIEF PRECISE COMPRESSION PEAK RECOVERY


BigHorn Ice Plus delivers cold therapy with pneumatic compression for faster recovery. Whether recovering from surgery, managing pain, or pushing athletic limits, its adjustable pressure system improves circulation while delivering consistent cooling. With precise temperature and timer controls, BigHorn Ice Plus enhances recovery, improves mobility, and keeps you moving.

HEAL FASTER WITH LESS PAIN AND MORE COMFORT


- ❏ Reduces Pain and Swelling
- ❏ Speeds Up Healing
- ❏ Comfortable, Flexible Wrap
- ❏ Whisper-Quiet Operation
- ❏ Pneumatic Compression Therapy
- ❏ Portable and Easy to Use
- ❏ Adjustable Temperature, Timer, and Compression




3 EASY STEPS TO TARGETED RELIEF



STEP 1
Fill the ice bucket with water and ice to the marked levels.



STEP 2
Connect the wrap and add the barrier.



STEP 3
Secure the wrap and press power to begin therapy.

General

Product Dimensions Ice Bucket: (±2mm) 235mm (L) x 235mm (W) x 280mm (H)

Cooling Duration ~2 hours per charge of battery, or directly plug into the wall outlet

Power Source Internal 5V 1A adapter

Volume 6.2L

Charging Time 4-5 hours

Adjustable Parameters 20-60 minutes of timer, 45-55°F of temperature, 0-60 mmHg of compression

Waterproof Grade IP20

Product Name	Part Number	Product Name	Part Number
Ice Bucket	(Tru- U102) *1 piece	Universal Pad	(Tru-PP102) *1 piece

Manual	(Tru-M102) *1 piece	Elastic Straps	(Tru-ES102) *3 pieces
--------	---------------------	----------------	-----------------------

AC Adapter	(Tru- AC102 Input: AC100-240V; Output: DC5.0V-1A) *1 piece	Non-Woven Fabric Barrier	(Tru-FB102) *1 piece
------------	--	--------------------------	----------------------



**TAKE CONTROL OF
YOUR RECOVERY**

GET STARTED